Table of Contents

Important Safety Instructions	03
Alternative Site Testing (AST)	06
Important Information about Your System	08
Intended Use.....	08
Test Principle.....	09
Introduction of Your System	10
System Contents.....	10
Meter Components.....	11
Meter Display.....	12
Test Strip.....	13
Meter Set-Up	13
Install Batteries.....	13
Set the Meter.....	15
Use the Alarm.....	22
Cancel the Alarm.....	22
Speaking Feature.....	23
Quality Control	24
Important Control Solution Test Information.....	25
Run a Control Solution Test.....	26
Blood Glucose Test	29
Prepare for the Test.....	29
Run a Blood Glucose Test.....	31

Expected Blood Glucose Values.....	35
Limitations.....	36
Comparison of Your Meter and Laboratory Results.....	37
Meter Memory.....	39
7-, 14-, 28-, 60- Day Average Results.....	39
Test Results.....	41
Your System Maintenance.....	44
Clean and Disinfection.....	44
Storage.....	53
Message and Troubleshooting.....	54
Specifications.....	61
Statements.....	62
Warranty.....	63
Appendix:	
Log Book	
Warranty card	

Important Safety Instructions

- The meter and lancing device are for single patient use. Do not share them with anyone including other family members! Do not use on multiple patients! *^{1,2}
- All parts of the kit are considered biohazardous and can potentially transmit infectious diseases, even after you have performed cleaning and disinfection procedures. *^{1,2} Details of the cleaning and disinfecting procedure could be found in the section Cleaning and Disinfection in this manual.
- Due to the potential biohazardous concern, you should wash hands thoroughly with soap and water after handling the meter, lancing device, or test strips.
- Severe dehydration and excessive water loss may cause false low results.
- If you believe you are suffering from severe dehydration, consult your healthcare professional immediately.
- If your blood glucose results are significantly lower or higher than your usual test results, and you do not have symptoms of illness, please repeat the test. If you have symptoms or continue getting results significantly higher or lower than your usual test results, please follow the treatment advice of your healthcare professional.
- If you are experiencing symptoms that are not consistent with your blood glucose test results and you have followed all instructions described in this Owner's Manual, please call your healthcare professional.

- Inaccurate results may occur in severely hypotensive individuals or patients in shock. Readings which are lower than actual values may occur for individuals experiencing a hypoxia state, or a hyperglycemic hyperosmolar state, with or without ketosis. Critically ill patients should not be tested with blood glucose meters.
- **Lipemic Effects:** Cholesterol up to 500 mg/dL or triglycerides up to 2000 mg/dL do not significantly affect the results. However, the glucose test may be affected beyond these levels.
- **Interferences:** When tested following NCCLS guidelines*³, bilirubin , creatinine, methyldopa ,galactose, maltoase ,xylose, ibuprofen, and salicylate at normal therapeutic concentrations do not significantly affect glucose results. However, these levels of the following interferents in blood may cause inaccurate test results:
 - Acetaminophen $\geq$ 12.5 mg/dL (Therapeutic level is 1.2-3.6 mg/dL)
 - Ascorbic acid $\geq$ 7.5 mg/dL (Therapeutic level is 0.4-2.1 mg/dL)
 - Dopamine $\geq$ 3 mg/dL (Therapeutic level is 0.04 mg/dL)
 - L-dopa $\geq$ 4 mg/dL (Therapeutic level is 0.02-0.3 mg/dL)
 - Tolbutamide $\geq$ 150 mg/dL (Therapeutic level is 3.6-7.2 mg/dL)
 - Uric acid $\geq$ 20 mg/dL (Normal level is 2.6-7.2 mg/dL)
 - Gentisic acid $\geq$ 25 mg/dL (Therapeutic level is 0.2-0.6 mg/dL)
 - Tolazamide $\geq$ 15 mg/dL (Therapeutic level is 2.0-2.5 mg/dL)


- Mannose $\geq$ 250 mg/dL (Therapeutic level is 1.15 mg/dL)
- Ibuprofen $\geq$ 50 mg/dL (Therapeutic level is 1.0-7.0 mg/dL)

References

1. FDA Public Health Notification: Use of Fingerstick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens: Initial Communication (2010)
<http://www.fda.gov/MedicalDevices/Safety/AlertsandNotices/ucm224025.htm>
2. CDC Clinical Reminder: Use of Fingerstick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens (2010)
<http://www.cdc.gov/injectionsafety/Fingerstick-DevicesBGM.html>
3. NCCLS, Vol. 25, No. 27, 2005, EP7-A2, Interference testing in Clinical Chemistry; Approved Guideline – Second Edition.
4. Tietz Textbook of Clinical Chemistry and Molecular Diagnosis (4th Ed.), Burtis, Ashwood & Bruns (Eds), Elsevier Saunders, 2005; 2301.
5. Bioavailability of Tolazamide from Tablets: Comparison of **In Vitro** and **In Vivo** results, P. G. Welling, R. B. Patel, U. R. Patel, W. R. Gillespie, W. A. Craig, K. S. Albert, 2006.

Alternative Site Testing (AST)

This system gives you to test on the palm, the forearm, the upper arm, the calf, or the thigh with the equivalent results to fingertip testing. These areas of your body are called alternative sites other than fingertip. Since nerve endings (receptors) are not so condensed at alternative sites, tests performed there will be less painful.

When you should use AST

Alternative site results may be different from fingertip results when glucose levels are changing rapidly.

Therefore, **use AST only when it has been more than two hours after any one of the following:**

- A meal.
- Diabetes medication.
- Exercise.

When testing blood glucose during or immediately after meals, medication, physical exercise, or any other event, **take blood sample from your finger only.**

Do NOT use blood samples from AST when

- You think your blood glucose is low.
- You have been diagnosed with hypoglycemic unawareness.
- The AST results do not agree with the way you feel.
- You are ill or stressed.
- It is within two hours of eating a meal, taking insulin or exercising.

AST results serve as references only. Consult your healthcare professional before performing alternative site testing (AST).**Before you get a drop of blood for AST**

Rub the puncture site about 20 seconds before penetration for increasing the accuracy.

Important Notices

- Alternative site testing can be used only during steady-state blood glucose conditions.
- Measurements from alternative site testing should never be used to calibrate a continuous glucose monitor (CGM) or entered into insulin dose calculators for insulin dosing recommendations.
- Alternative site testing should not be used for insulin dose calculations.

Important Information about Your System

Intended Use

ADVOCATE® Redi-Code+ BMB-EA001S Blood Glucose Monitoring System is intended for use outside the body (in vitro diagnostic use) at home. It is used for quantitative measurement of glucose level in fresh capillary whole blood samples (from the finger, the palm, the forearm, the upper arm, the calf, or the thigh). The alternative site testing can be only used during steady-state blood glucose monitoring. The ADVOCATE® Redi-Code+ BMB-EA001S Blood Glucose Monitoring System is intended for use by a single person and should not be shared. In addition, it is intended for use at home as an aid in monitoring the effectiveness of diabetes control program. It should not be used for the diagnosis of diabetes, or for the testing of neonates.

The ADVOCATE® Redi-Code+ BMB-EA001S Blood Glucose Monitoring Systems consists of the ADVOCATE® Redi-Code+ BMB-EA001S Blood Glucose meter and the ADVOCATE® Redi-Code+ BMB-BA006A Blood Glucose test strips. The ADVOCATE® Redi-Code+ BMB-EA001S meter are used only with ADVOCATE® Redi-Code+ BMB-BA006A Blood Glucose test strips to quantitatively measure glucose in fresh capillary whole blood samples drawn from finger tips, the palm, the forearm, the upper arm, the calf, and the thigh. **The meter also includes some speaking functions but has not been validated for use by the visually impaired.**

ADVOCATE® Redi-Code+ control solutions are for use with the ADVOCATE® Redi-Code+ BMB-EA001S Blood Glucose Monitoring System as a quality control check to verify the accuracy of blood glucose test results.

Test Principle

ADVOCATE® Redi-Code+ BMB-EA001S Blood Glucose Test is based on the quantitative measurement of glucose level in whole blood. The blood sample is drawn into the tip of the test strip for testing. Glucose in the sample mixes with special chemicals on the test strip and a small electrical current is produced. The current is measured by the meter and the corresponding blood glucose level is displayed. The strength of the current produced by the reaction depends positively on the amount of glucose in the blood sample.

Cautions:

- Only use ADVOCATE® Redi-Code+ BMB-BA006A Test Strips with ADVOCATE® Redi-Code+ BMB-EA001S Meter. Other test strips may produce inaccurate results.
- Perform all quality control checks recommended in this User's Manual.
- As with all small parts in this system, these parts should be kept away from children who still put things in their mouths. If they are swallowed, promptly see a healthcare professional.
- Read all instructions thoroughly and practice the test before using the product to test your blood glucose. Do all quality control checks as directed and consult with a diabetes healthcare professional.

Introduction of Your System

System Contents

KIT CONTENTS		
Meter	Owner's Manual	Quick Reference Guide
		

English

Notes:

- If this system has been opened prior to use, please return it to the place of purchase.
- For information on how to obtain accessories including test strips, lancets and control solutions. Please contact the place of purchase for availability.

Meter Components

Test Strip Port

Insert test strip. The meter turns on when you insert the test strip.

Display


Display results and other important information.

Main Button

Press to enter memory mode to review stored results or adjust setting.



Set Button


Press to enter setting mode and change the date, time, alarm settings and memory deletion.

Battery Cover

Remove cover and replace battery.

Test Strip

Confirmation Window

Blood is drawn in and reacts with enzyme.

Contact Bars



Insert this area to test strip port.

Meter Set-Up


Install Batteries

Make sure that the meter is turned off before replacing the batteries.

1. Press the buckle on battery cover and lift up to remove cover.
2. Remove the old battery and replace with two 1.5V AAA alkaline batteries.



3. Close the battery cover.


Notes:

- Batteries should be replaced when the power of your meter is getting low by displaying battery symbol or the meter briefly shows “E-2” with battery symbol.

- After replace new batteries, the display check will appear briefly. That means the batteries are installed successfully.

- Replacing the batteries does not affect the meter's memory; however, the meter's setting (date, time and alarms) may need to be re-set.
- Batteries might leak chemicals if not used for a long time. Remove the batteries if you are not going to use this system for an extended period (i.e., 3 months or more).

Set the Meter

Your meter allows you to set the year, month, date, time, alarm as well as delete memory.

- **Enter Setting Mode**

Start with the meter off. Press Set button (the back of the meter).

- **Exit Setting Mode**


Any time you want to exit setting mode, press and hold Set button for 2 seconds until the meter turns off.

- **Move Setting Faster**

Hold M button down to scroll forward quickly until the number you desire appears.


1. Set the Year

Press Set button to enter Setting mode. The Year will be flashing. Press and release M button until the correct year appears. Press Set button to move to Month setting.

2. Set the Month

With the month flashing, press and release M button until the correct month appears. Press Set button to move to Day setting.


3. Set the Day

With the day flashing, press and release M button until the correct day appears. Press Set button to move to Hour setting.


4. Set the Hour

With the hour flashing, press and release M button until the correct hour appears. Press Set button to move to Minutes setting.


5. Set the Minutes

With the minutes flashing, press and release M button until the correct minute appears. Press Set button to move to Alarm setting.

6. Set the Alarm

The meter provides six alarms (1-6). You may set up any or all of them. With the first alarm flashing (the default is OFF), press M button to switch "OFF" to "ON" and go to **step 7-1** to see how to set the alarm's time. Or with "OFF" displaying, keep pressing Set button to move to Language setting.

7-1. Enter the Alarm Setting

01 (Alarm 1) and 📶 with flashing “On” displays, press Set button to enter the hour setting.

7-2. Set the Alarm's Hour

With the hour flashing, press and release M button until the correct hour appears. Press Set button to move to Minutes setting.


7-3. Set the Alarm's Minutes

With the minutes flashing, press and release M button to select 15 minutes, 30 minutes, or 45 minutes. Press Set button to finalize Alarm 1 setting and then move to Alarm 2 setting. If you want to leave alarm setting, when “OFF” displays with alarm marker, press Set button to move to Language setting.


7-4. Set Other Alarms (2-6)

Repeat steps 7-1, 7-2 and 7-3 to set Alarm 2, Alarm 3, Alarm 4, Alarm 5 and Alarm 6. After finalizing Alarm 6, press Set button to move to Language setting.


8. Select the Speaking Language

The default language is English (indicated by flashing "E"). Press M button to toggle the language setting to Spanish (indicated by flashing "S") and press M button again to toggle it back. Press Set button to move to Volume setting.


9. Set Speaking Volume

Press M button to adjust the speaking volume(0 for the silent mode and 6 for the largest speaking volume). Press Set button to move to Memory Deletion.


10. Delete Memory

With "M" symbol and flashing "dEL" appears, if you **DO NOT** want to delete memory, press Set button to skip this step. If you want to delete ALL memory, press M button for 2 seconds to clear ALL memory. With "CLr" displaying, that means ALL memory is deleted.

After pressing Set button or deleting memory, set-up is completed. The meter will show "OFF" and then turn off automatically.


OFF

Notes:


- The date, time and alarm can **ONLY be changed** in the setting mode. Therefore, when performing a glucose testing, those parameters are not possible to be changed.
- The meter displays 7-, 14-, 28-, and 60-day averages which you can access from the meter memory. These averages are calculated from results obtained during the 7-, 14-, 28-, and 60-days preceding the current date and time settings. When the date and time are changed, the 7-, 14-, 28-, and 60-days averages will be re-calculated and updated automatically based on the new date.
- The date and time might need to be adjusted due to daylight saving.
- The date, time, language and voice volume setting might need reset after you replace the battery.
- If no button is pressed for 2 minutes, the meter will turn off automatically.

Use the Alarm

When the alarm sounds, the meter will automatically turn on. The alarm will keep sounding for 30 seconds. You can perform a test or press M button to turn off the meter.

**Perform a Test**

Insert a test strip into the test strip port and follow the Blood Glucose Test instructions (page 29).

**Turn off the Alarm**

Press M button to turn off the meter or leave the alarm sounding for 30 seconds and then the meter will turn off.

Cancel the Alarm

If you want to cancel the alarm, follow manual's "Set the Meter" instructions.

Speaking Feature

The ADVOCATE® Redi-Code+ Blood Glucose Meter has designed to be with speaking function. With the speaking function on, you will be greeted by the meter every time you turn on the meter. The speaking function, determined by which function you are currently operating, will guide you through the next action, or state your current mode and record. These features are to aid you in operating the meter, but it has not been validated for use by the visually impaired. If you wish to turn the function off, please set the volume to 0 by referring to “Set Speaking Volume” subsection in “Meter Set Up”. If you wish to change the language from preset English to Spanish, please to refer to “Select the Speaking Language” subsection in “Meter Set Up”.

* The speaking function is for additional reference only, please refer to the meter screen for you blood glucose measuring result.

Quality Control

ADVOCATE® Redi-Code+ Level 1/Low and Level 3/High Control Solutions are able to check the meter and the test strips are working together. It is very important that you do this simple check routinely to make sure your results are accurate.

When to check your system with control solutions

- Before using the meter for this first time.
- When you begin using a new vial of test strips.
- Whenever you suspect that the meter or test strips are not working properly.
- After cleaning and disinfection of the meter.
- At least once a week.
- To check your technique.
- When your blood glucose test results are not consistent with how you feel, or when you think your results are wrong.
- When your test strips are exposed to extreme environmental conditions (below 39.2°F, or above 104°F, or humidity levels above 85%).
- If you drop the meter.

Important Control Solution Test Information


The control solution range printed on the test strip vial is for ADVOCATE® Redi-Code+ Control Solution only. It is used to evaluate meter and test strip performance. It is not the recommended range for your blood glucose level.

- Use only ADVOCATE® Redi-Code+ control solutions. Check the expiration date on the control solution vial. Do not use if expired.
- Control solution, meter, and test strips should come to room temperature between 68°F and 77°F before testing.
- Shake the vial before use. Discard the first drop of control solution. After squeezing, wipe off the dispenser tip to avoid contamination. The above tips ensure you get a good sample and an accurate result.
- Due to the potential biohazard concern, you should wash hands thoroughly with soap and water after handling the meter, lancing device, or test strips.
- Use for only 90 days from the date that you first open it. Record the discard date (date opened plus 90 days) on the control solution vial.
- Discard after 90 days.
- Store the control solution tightly closed at temperatures between 35°F - 86°F.
- Do not freeze.

Run a Control Solution Test


1. Insert a Test Strip

Fully insert a test strip to the meter with the contact bars end first. Push it in until it will go no further. The meter will turn on and the display check will appear briefly.


2. Press M Button for CTL mode

When the flashing apply blood symbol appears on the display, press M button thrice until a control bottle with "CTL" appears. **Control results marked with a control bottle and "CTL" will NOT be stored in the meter's memory.**

3. Obtain Control Solution


Shake the control solution vial well. Remove the cap and squeeze the vial to discard the first drop. Then wipe off the dispenser tip to prevent contamination. Squeeze the vial again to produce another drop and **place this drop on a clean surface.**

4. Apply Control Solution


Hold the meter to touch the drop until the meter beeps and begins to count down.

To avoid contaminating the control solution with the content of the test strip, DO NOT directly apply control solution into a strip.

5. Read the Control Solution Result

After 5 seconds, the control solution test result appears. Compare the result with the range printed on the test strip vial. The result should fall within this range. Remove the used strip, the meter will automatically turn off.

English

Cautions:

- The control solution range printed on the test strip vial is for ADVOCATE® Redi-Code+ Control Solution only. It is not a recommended range for your blood glucose test results.
- Every time you perform a control solution test, you must enter into the “CTL” mode so that the result will not be stored in the meter memory. Failure to do so will confuse the blood glucose result with the control solution result in meter's memory.

If the result is not in the acceptable range, it may be caused by:

- Performing the test improperly.
- Not shaking the control solution vial well enough.
- Using expired or contaminated control solution and test strips.
- Using the first drop of control solution.
- Using a damaged meter.
- Control solution is stored at environment out of the specified temperature range, which is from 35°F to 86°F.

If you continue to get control solution test results that fall outside the range, the system may be not working properly. **Do NOT** use the system to test your blood. **Contact your Customer Support at 1-877-394-1655 for help.**

Blood Glucose Test

Make sure you have all items needed to test before you perform your blood glucose test.

- Blood Glucose Meter
- Test Strip
- Lancing Device
- Sterile Lancet

Warnings:


To reduce the chance of infection

- Make sure to wash and dry the puncture site before starting.
- Never share a lancet or the lancing device with anyone to avoid the risk of blood borne pathogen transmission.
- Always use a new, sterile lancet. Lancets are for single use only.
- Keep the meter, the lancets and the lancing device clean before using.
- The user should wash hands thoroughly with soap and water after handling the meter, lancing device, or test strips.
- If the meter is being operated by a second person who is providing testing assistance to the user, the meter and lancing device should be decontaminated prior to use by the second person. Please refer to the section "Cleaning and Disinfection" for details in this manual.

Prepare for the Test

1. Insert a Lancet

Remove the cap and then insert a lancet into the holder and push in firmly.

2. Remove the Protective Disk

Twist the protective disk until it separates from the lancet and leave the cap aside for later use.

3. Replace the Lancing Device Cap

Replace the lancing device cap by twisting it back on until it is snug.


4. Set the Lancet Depth

Twist the adjustable tip in either direction until the number lines up with the arrow. Start at Level 3 or 4, for tougher skin, twist to a higher number.

5. Cock the Lancing Device

Slide the ejection/cocking control back until it clicks. If it does not click, it may have been cocked when the lancet was inserted.

The lancing device is now ready for use.

Run a Blood Glucose Test

1. Insert a Test Strip

Fully insert a test strip to the meter with the contact bars end first. Push it in until it will go no further. The meter will turn on and the display check will appear briefly.

2. Choose Pre-Meal/Post-Meal/General Mode

You can choose three modes by pressing M button, if you miss the desired mode, you can continue pressing M button to go back.

Pre-Meal Mode


When the flashing apply blood symbol appears on the display, pre-meal marker (☞) displays firstly. Go to Step 3 for obtaining a drop of blood.

Post-Meal Mode

When the flashing apply blood symbol appears on the display, press M button once and postmeal marker (☛) displays. Go to Step 3 for obtaining a drop of blood.

General Mode

When the flashing apply blood symbol appears on the display, press M button twice and general mode appears (no marker displays). Go to Step 3 for obtaining a drop of blood.

3. Obtain a Drop of Blood


Before puncturing, clean the puncture site with 70% alcohol cotton and let it air-dry.

Fingertip

Hold the lancing device firmly against the side of your fingertip and press the release button. After hearing a click, it indicates that the puncture is complete.

Sampling for an alternative site (AST) other than fingertip

Press and hold the lancing device against your palm, forearm, upper arm, calf or thigh for a few seconds, then press the release button. Keep holding the lancing device and cap against your skin until a round drop of blood appears.

Warnings:

- The risk of transmitting of disease such as viral hepatitis from shared use of lancing device has been increasingly reported by the FDA and CDC ^{*1,2}. To prevent the transmission of disease and infection, the lancing device is intended only for a single user and should not be shared.
- Squeeze the punctured area gently to obtain blood and be careful NOT to smear the blood sample.
- Choose a different spot each time you test. Repeated punctures in the same spot may cause soreness and calluses.
- Before you decide to sample for alternative sites, please refer to "Alternative Site Testing" for limitations or consult your healthcare professional.
- It is recommended that you discard the first drop of blood as it might contain tissue fluid and serum, which may affect the test result.

4. Apply Blood to the Test Strip

Apply the tiny blood sample to the test strip tip. The minimum blood sample size is 1.1 μ l (●). Blood will be drawn into the strip. Keep holding the drop of blood to the hole of the test strip until the confirmation window is full and the meter begins to count down.

Please note that if the confirmation window is not filled with blood, the test result may be inaccurate.



5. Get Your Result

After 5 seconds, your blood glucose test result appears. This result will be automatically stored in the meter's memory.

6. Dispose the Used Lancet

Remove the lancet carefully to prevent injury.

Take the lancet out from the lancing device. Place the disk on a hard surface and push the exposed tip into the disk cover.

English

Warnings:

- The user should wash hands thoroughly with soap and water after handling the meter, lancing device, or test strips.
- If the meter is being operated by a second person who is providing testing assistance to the user, the meter and lancing device should be decontaminated prior to use by the second person. Please refer to the section “Cleaning and Disinfection” for details in this manual.
- If no action is done for 2 minutes, the meter will automatically turn off. You must remove and reinsert the test strip to restart the test procedure.
- The blood should completely fill the confirmation window before the meter begins to count down. If you find that the confirmation window is not filled with blood when the meter is counting, **NEVER** try to add more blood to the test strip. **Discard the test strip and retest with a new one.**
- Any trouble filling the confirmation window, contact your Customer Support at 1-866-373-2824 for help.
- The used lancet and the used test strip may be potential biohazard. Please discard it carefully according to your healthcare provider's instructions.

Expected Blood Glucose Values

Many factors will affect blood glucose values such as food intake, medication dosages, health, stress, or exercise. The expected blood glucose values appropriate for you should be consulted and determined by your healthcare professional. Blood glucose levels for people with diabetes are expected to be controlled in the following ranges*⁶ :

<i>Time</i>	<i>Range (mg/dL)</i>
Fasting and before meals	70-130
1- 2 hours after meals	Less than 180

Note: The normal blood glucose level for people without diabetes is 126 mg/dL.*⁵

References

6. American Diabetes Association (2010). Standards of Medical Care in Diabetes. Diabetes Care, 33 (Supplement 1): S11–S61.

Limitations

- a. The test strips are used for fresh capillary whole blood samples.
- b. DO NOT use neonate blood sample.
- c. Not to be used for diagnosis or screening of diabetes.
- d. Alternative site testing with this system can be used only during steady-state blood glucose conditions.
- e. Measurements from alternative site testing should never be used to calibrate a continuous glucose monitor (CGM) or entered into insulin dose calculators for insulin dosing recommendations.
- f. Extreme humidity may affect the results. A relative humidity greater than 85% may cause incorrect results.
- g. The system should be used at a temperature between 50°F - 104°F (10°C - 40°C). Outside this range, the system halts.
- h. DO NOT reuse the test strips. The test strips are for single use only.
- i. Red blood cell count (Hematocrit); Hematocrit levels between 20% and 60% will not significantly affect the results, Hematocrit levels below 20% may cause higher results, Hematocrit levels above 60% may cause lower results. If you do not know your hematocrit level, please consult with your healthcare professional.
- j. Altitude up to 10,745 feet above sea level has no effect on readings.

Comparison of Your Meter and Laboratory Results

Test results with this meter are plasma-calibrated. Since the laboratory tests provide plasma-calibrated results, your meter may differ from your laboratory result due to normal variation. Meter results can be affected by factors and conditions that do not affect laboratory results in the same way. To obtain an accurate comparison between meter and laboratory results, follow a few basic guidelines below.

Before going to the lab

- Perform a control solution test to make sure that the meter is working properly.
- Fast for at least eight hours before you perform your blood test.
- Take your meter with you to the lab.

While at the lab

- Only use fresh capillary blood obtained from the fingertip.
- Collect your blood into heparin-containing test tubes and do **NOT** use other anticoagulants or preservatives.
- Conduct your meter test within 15 minutes of the lab test.

You may still have a variation from the result because blood glucose levels can change in some specific situations that could cause a difference of more than $\pm 20\%$.

- You have eaten recently. The blood glucose level from a finger tip can be up to 70 mg/dL higher than blood drawn from a vein (venous sample) used for a lab test^{*7}.
- If your hematocrit (the amount of red blood cells in the blood) is high or low.
- You are suffering from severe dehydration (the loss of body fluid).

You may see the package insert of test strip for accuracy and precision data and for important information on limitations.

References

7. Sacks, D.B.: "Carbohydrates." Burtis, C.A., and Ashwood, E.R.(ed.), Tietz Textbook of Clinical Chemistry. Philadelphia: W.B. Saunders Company (1994), 959.

Meter Memory

The meter stores the last 400 blood glucose test results along with the date and time. It also provides you with 7-, 14-, 28- and 60-day averages of your pre-meal, post-meal and general blood glucose test results, respectively. You can begin to review the results stored in meter's memory by following the next steps.

Enter Memory Mode

Start with the meter off. Press M button until a test strip marker appears.

Exit Memory Mode

Any time you want to exit memory mode, press and hold M button for 2 seconds until the meter turns off.

7-, 14-, 28-, 60- Day Average Results

1. Enter Memory Mode

When the meter is off, press and release M button.

The meter will show flashing strip symbol, it means you can test your blood glucose right now. If you do not want to test a blood, press M button to review 7-day average result for pre-meal mode.


2. 7-Day Average Pre-Meal Result

The 7-day average pre-meal result appears after press and release M button.

07 DAY AVG = your 7-day average

30 = indicates the number of results included in your average

 = pre-meal marker


3. 7-Day Average Post-Meal Result

Press and release M button, then 7-day average post-meal result appears.

07 DAY AVG = your 7-day average

30 = indicates the number of results included in your average

— = post-meal marker


4. 7-Day Average General Result

Press and release M button, then 7-day average general result appears.


07 DAY AVG = your 7-day average

30 = indicates the number of results included in your average


5. 14-Day Average Results for Pre-Meal, Post-Meal and General Mode

Press M button again, 14-day average pre-meal result appears. Like Step 3 to 4, 14-day average post-meal result and general result will appear in order.


6. 28-Day Average Results for Pre-Meal, Post-Meal and General Mode

After 14-day average general result, press M button, 28-day average pre-meal result appears. Like Step 3 to 4, 28-day average post-meal result and general result will appear in order.

7. 60-Day Average Results for Pre-Meal, Post-Meal and General Mode

After 28-day average general result, press M button, 60-day average pre-meal result appears. Like Step 3 to 4, 60-day average post-meal result and general result will appear in order.


Continue to press the M button to view stored test results. The most recent result will display first.

Test Results


1. View Test Results

After viewing 60-day average result, press M button again. The most recent result with date and time will be shown on the display.

2. View Second Most Recent Result

Press the M button one by one to recall the test results stored in the meter consecutively.

3. View the Last Results



After viewing the last result, press M button again and the meter will be turned off.

Notes:

- The control solution results are NOT stored in the memory. The list of past results and the result average are for blood glucose results only.
- The stored results marked with a pre-meal marker (⏏)/post-meal marker (⏏) depend on what you chose in your test procedure. If you chose general mode in your test procedure, the stored results without markers will display.

- If no button action is detected for 2 minutes, the meter will show “OFF” and turn off automatically.
- If no results are stored in the meter's memory, “---” will be shown on the LCD.

Your System Maintenance

Cleaning and Disinfection

The risk of transmission of disease such as viral hepatitis from shared use of the meter and lancing device has been increasingly reported by the FDA and CDC^{*1,2}. The cleaning procedure is to thoroughly remove all debris and bioburden from surfaces and the disinfection procedure is to destroy or prevent growth of organisms capable of infection, **the**

ADVOCATE® Redi-Code+ BMB-EA001 Blood Glucose Meter and Lancing Device should be cleaned and disinfected when the following situations happen.

- **Whenever the meter and lancing device are visibly dirty.**
- **At least once per week.**
- **Before the meter and lancing device are being operated by a second person who is providing testing assistance to the user.**


Make sure you have all items before you clean and disinfect the meter and lancing device.

- A commercially available disinfectant towelette, CaviWipesXL (manufactured by Metrex with EPA registration number 46781-8). Please contact Metrex Research Corp. at 1-800-841-1428 for availability. The CaviwipesXL is also available via on line purchase from website such as <http://www.doctordepotonline.com/>


- The CaviWipesXL has been shown to be safe for use with ADVOCATE® Redi-Code+ BMB-EA001S Blood Glucose meter and lancing device, but any other disinfectant product with the EPA registration number may be used on this device.

* Please follow recommendations in the CaviWipesXL (EPA Reg. No. 46781-8) labeling for the safe use of this product.

Area to be cleaned and disinfected for the meter

Area to be cleaned and disinfected for the lancing device


Do	Do Not
• Make sure the meter is turned off during cleaning and disinfection. • Keep the test strip vial tightly closed when performing the cleaning and disinfection procedures because the fumes from the disinfectant may affect the performance of the strip. • After cleaning or disinfection, please perform the physical appearance and performance check of devices.	• Get any moisture in the test strip slot. • Spray any cleaning solution directly onto the meter. • Put the meter under water or liquid. • Pour liquid into the meter.

1. Cleaning Process

1-1. Put on disposable protective gloves (based on the “Instruction of Use” of CaviwipesXL) and take disinfectant CaviWipesXL towelette from the container.

1-2. Wipe down the body of the meter and lancing device, until the debris or the bioburden is removed.
After air-drying, go to Step 2, Disinfection Process for the meter and lancing device.


1-3. Discard the used disinfectant wipe in trash.

2. Disinfection Process

* A pre-cleaning step (using the CaviWipesXL) is required prior to each disinfection step.

2-1. Take a new disinfectant CaviWipesXL from the container.

2-2. Wipe down the body of the meter and lancing device enough to thoroughly wet the surface. Repeat this step for six times (3 passes horizontal and 3 passes vertical as shown in the figure) to complete the disinfection procedure. Allow the devices to remain wet for a 2-minute contact time. Similarly clean and disinfect the lancing device, including the clear cap used for alternative site testing. Do disinfection at least once per week.

2-3. Discard the used disinfectant wipe in trash.

2-4. The user should wash hands thoroughly with soap and water after handling the meter, lancing device, or test strips.

After each cleaning or disinfecting process, you should perform the physical appearance check and performance check of the meter and lancing device for items shown in the tables below.

Physical appearance check-after each cleaning and disinfection for the meter

Check item	Acceptance Standard
Losing display transparency	Accept the result if none of the items listed on the left column appears after each cleaning or disinfection process
Losing legibility of any marking	
Polymer crazing	
Cracking	
Swelling	
Dissolving	
Softening	
Becoming brittle	
Deformation	
Cloudy	

Performance check-after each cleaning or disinfection for the meter

Steps	Check items	Accept result
1. Insert test strip	Is the meter powered on?	Accept result if the meter is powered on
2. Press M Button for CTL mode	Is "CTL" showing on the display?	Accept result if the "CTL" is showing on the display
3. Apply control solution to obtain the control solution test result	Is the control test result within the control solution range specified on the strip vial?	Accept result if the control test result is within the control solution range specified on the strip vial

* Stop using the meter if signs of deterioration are noted.

Physical appearance check- after each cleaning or disinfection for the lancing device, including lancing device cap

Check item	Acceptance Standard
Polymer crazing	Accept the result if none of the items listed on the left column appears after each cleaning or disinfection process
Cracking	
Swelling	
Dissolving	
Softening	
Becoming brittle	
Deformation	
Cloudy	

Performance check-after each cleaning or disinfection for the lancing device cap


Steps	Check items	Accept result
1. Lancing device cap	Lancing device cap. Fit perfectly with the main body of the device?	Accept result if both caps fit perfectly with the main body of the device
2. Cock the device and press on the trigger	Does the device work after pressing the trigger?	Accept result if the device works
3. Install a lancet to the lancet holder and push in firmly	Could the lancet holder hold the lancet firmly?	Accept result if the holder could hold the lancet firmly

* Stop using the lancing device if signs of deterioration are noted.

In case the result is not acceptable, please contact Customer Support service, with the contact information listed on the last page of the manual, for replacement with a new device.

Note:

- **Special attention should be given during the cleaning and disinfecting procedures- DO NOT squeeze the towel, which might cause meter malfunction by soaking the test strip port with excess liquid disinfectant.**

Virucide cannot drop into the strip port.



- You should wash hands thoroughly with soap and water after handling the meter, lancing device, or test strips.
- Any trouble performing cleaning and disinfection procedure on the meter and lancing device, **please contact your Customer Support at 1-866-373-2824 for technical assistance.**
- Stop using the lancing device if signs of deterioration are noted.
- Based on the robustness study performed by the manufacturer, the meter has passed all the appearance and performance checks after a total of 3,000 cleaning and disinfection cycles (over 5 years), and the lancing device has passed all the appearance and performance checks after a total of 5,000 cleaning and disinfection cycles (over 5 years). The meter must be replaced with a new meter at the end of 3,000 cleaning and disinfection cycles or 5 years, whichever comes first. The lancing device must be replaced with a new lancing device at the end of 5,000 cleaning and disinfection cycles or 5 years, whichever comes first.



Storage




- Always store or transport your meter, test strips, and control solution in the carrying case after each use.
- Store your meter, test strips, and control solution in a cool and dry place. Keep away from direct sunlight and heat.
- Avoid dropping and strong impact on your system.
- After removing a test strip from the vial or using control solution, immediately replace the vial cap and close it tightly to avoid contamination or damage.
- Do not use test strips and control solution after expiration date. Write the discard date on their vials after first opening.
- Keep the strip vial away from children since the cap and the test strip may be a choking hazard. If swallowed, promptly see a healthcare professional for help.




Messages and Troubleshooting

This system displays messages when there are problems with the test strip, with the meter, or when your blood glucose levels are lower than 20 mg/dL or higher than 600 mg/dL. These messages help to identify certain problems but do not appear in all cases when a problem has occurred. Improper use may cause an inaccurate result without producing an error message. In the event of a problem, refer to the information under “What to do”.

Message	What it means	What to do
	Your blood glucose result is lower than 20 mg/dL (below measurement limit).	You should seek immediate medical assistance.
	Your blood glucose result is higher than 600 mg/dL (above measurement limit).	You should seek immediate medical assistance.

Message	What it means	What to do
	Your blood glucose result is equal to or higher than 240 mg/dL, KETONE? appears with your result.	Please seek immediate medical assistance.
	✓ The test strip cannot be identified. ✓ A previously used test strip may have been inserted.	Use a new test strip to perform test strip. Be sure you use ADVOCATE® Redi-Code+ BMB-BA006A test strips.
	Insufficient blood sample or the test strip is not completely filled with blood sample.	Remove the test strip and retest using a new test strip. Be sure you apply sufficient blood sample to fill test strip.

Message	What it means	What to do
 The image shows a digital display with the text 'E-2' in large characters. Below the text is a small icon of a battery with a cross through it, indicating a dead battery.	Batteries are dead.	Replace the batteries immediately.
 The image shows a digital display with the text 'E-1' in large characters. Above the text is a small 'M' in a box. To the left of the text is a speaker icon. To the right of the text is a small battery icon. Below the text is the text 'mg/dL'. The background of the display shows a date and time: '03-10 06:38 AM'.	When battery symbol displays with the test result, battery power level is low.	Replace the batteries as soon as possible.
 The image shows a digital display with the text 'E-3' in large characters.	Operating temperature is out of range.	Move to an area between 50°F and 104°F, let the meter adjust to this temperature for 20 minutes before performing a test.

Message	What it means	What to do
	No results are stored in the meter's memory.	No action required.
	No average results are stored in the meter's memory.	No action required.
	The alarm sounds to remind you of performing a test.	Insert a new test strip to perform a test or press M button to turn off the meter.

Message	What it means	What to do
Meter does not display a message after inserting a test strip	Possible dead battery.	Replace the battery.
	Battery incorrectly installed or missing.	Check that the battery is correctly installed.
	Test strip inserted upside down or incompletely.	Insert the test strip correctly with the contact bars end.
	Defective meter.	Please call your Customer Support for help.

Message	What it means	What to do
Test does not start after applying the sample	Insufficient blood sample.	Repeat the test using a new test strip with larger volume of blood sample.
	Defective test strip.	Repeat the test with a new test strip.
	Sample applied after automatically shutoff (2 minutes after last user action).	Repeat the test with a new test strip. Apply sample only when blood drop symbol appears on the display.
	Defective meter.	Please call your Customer Support for help.

Message	What it means	What to do
Control solution test result is out of range	Error in performing the test.	Read the instruction thoroughly and repeat the test again.
	Do not shake the control solution vial very well.	Shake the control solution vigorously and repeat the test again.
	Expired or contaminated control solution.	Check the expiration date or the discarded date of the control solution.
	Control solution that is too warm or too cold.	Control solution, meter, and test strips should come to room temperature (68-77°F) before testing.
	Test strip deterioration.	Repeat the test with a new test strip.
	Defective meter.	Please call your Customer Support for help.
	Control solution is stored at environment out of specified temperature range, which is from 35°F to 86°F.	Use a new bottle of control solution that has been stored at specified temperature range.

Never try to disassemble the meter in any circumstances. If you encounter any error messages not listed above or if you have followed the actions recommended but the problem keeps unsolved, please call your Customer Support at 1-877-394-1655 for support.

Specifications

Model No.	BMB-EA001S
Dimension	64mmx95mmx29mm
Weight	53g without batteries
Power Source	Two 1.5V AAA alkaline batteries
Battery Life	Over 500 tests
Auto Turn-off	2 minutes after idle.
Test Sample	Capillary whole blood
Sample Size	1.1 microliter
Test Time	5 seconds
Measurement Range	20 to 600 mg/dL
Meter/Strip Storage Condition	39.2°F to 104°F, / 85% R.H. (noncondensing)
Operating Condition	50°F to 104°F, 15% / 85% R.H. (noncondensing)
Memory	400 blood glucose tests with date and time
Control Solution Storage Condition	35°F to 86°F

The device has been certified to meet the electrical and safety requirements of: IEC 61010-1, EN 61010-1, IEC 61010-2-101, EN 61010-2-101, IEC 61326-1, EN 61326-1, IEC 61326-2-6, EN 61326-2-6, IEC 60601-1-2, EN 60601-1-2

Statements

- The equipment complies with IEC 61326-1:2005. Emissions are low and not likely to cause interference in nearby electronic equipment.
- The equipment is tested for immunity to electrostatic discharge as specified in IEC 61000-4-2, at a level of 8kV for Air Discharge and 6kV for Contact Discharge.
- The equipment is tested for immunity to radio frequency interference at the frequency range of 80 MHz to 2.7GHz

Warranty

Diabetic Supply of Suncoast, Inc. warrants to the original purchaser that this instrument will be free from defects in materials and workmanship for a period on ONE (1) year from the date of purchase (except as noted below). The warranty is not transferable.

This warranty is subject to the following exceptions and limitations

1. This warranty is valid only if the Warranty Registration Card is completed with date of purchase and serial number and sent to Diabetic Supply of Suncoast's office within 14 days commencing from the date of purchase.
2. This warranty is limited to replacement due to defects in parts and workmanship.
3. This warranty does not apply to the performance of this product that has been damaged due to abuse, accidents, alteration, misuse, neglect, maintenance by someone other than Diabetic Supply of Suncoast, or failure to operate the instrument in accordance with instructions. Further, Diabetic Supply of Suncoast assumes no liability for malfunction or damage caused by the use of other manufactures redistributers other than Diabetic Supply of Suncoast.
4. Diabetic Supply of Suncoast reserves the right to make changes in design of this instrument without obligation to incorporate such changes into previously manufactured instruments.

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Referencias

1. FDA Public Health Notification: Use of Fingerstick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens: Initial Communication (2010)
<http://www.fda.gov/MedicalDevices/Safety/AlertsandNotices/ucm224025.htm>
2. CDC Clinical Reminder: Use of Fingerstick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens (2010)
<http://www.cdc.gov/injectionsafety/Fingerstick-DevicesBGM.html>
3. NCCLS, Vol. 25, No. 27, 2005, EP7-A2, Interference testing in clinical chemistry; approved guideline – second edition.
4. Tietz Textbook of Clinical Chemistry and Molecular Diagnosis (4th Ed.), Burtis, Ashwood & Bruns (Eds), Elsevier Saunders, 2005; 2301.
5. Bioavailability of tolazamide from tablets: Comparison of in vitro and in vivo results, P. G. Welling, R. B. Patel, U. R. Patel, W. R. Gillespie, W. A. Craig, K. S. Albert, 2006.

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




Mensaje	Qué significa?	Que hacer?
 The image shows a digital display with the text "E-2" in large characters. Below the text is a small battery icon with a single bar, indicating a low battery level.	Las baterías están agotadas.	Reemplace de inmediato las baterías.
 The image shows a digital display with a small "M" icon in the top left corner. The display shows the date and time "03-10 06:38 AM" and a large number "137" below it. Underneath the number is the unit "mg/dL". At the bottom left is a speaker icon, and at the bottom right is a battery icon with a single bar.	Cuando el símbolo de la batería se muestra junto con el resultado, el nivel de las baterías esta bajo.	Reemplace las baterías lo antes posible.
 The image shows a digital display with the text "E-3" in large characters.	La temperatura de operación está fuera del rango.	Muévase a una área entre 50°F and 104°F, permita a que el medidor se ajuste a esta temperatura durante 20 minutos antes de realizar una prueba.

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PHONE

My Self-Check Diary

	BREAKFAST				LUNCH			
WEEK	Pre-Meal 	Insulin Medication	Post-Meal 	Pre-Meal 	Insulin Medication	Post-Meal 		
MON								
TUE								
WED								
TUR								
FRI								
SAT								
SUN								
AVG								

My Self-Check Diary

	BREAKFAST			LUNCH		
WEEK	Pre-Meal 	Insulin Medication	Post-Meal 	Pre-Meal 	Insulin Medication	Post-Meal 
MON						
TUE						
WED						
TUR						
FRI						
SAT						
SUN						
AVG						

My Self-Check Diary


	BREAKFAST			LUNCH		
WEEK	Pre-Meal 	Insulin Medication	Post-Meal 	Pre-Meal 	Insulin Medication	Post-Meal 
MON						
TUE						
WED						
TUR						
FRI						
SAT						
SUN						
AVG						

My Self-Check Diary

	BREAKFAST				LUNCH			
WEEK	Pre-Meal 	Insulin Medication	Post-Meal 	Pre-Meal 	Insulin Medication	Post-Meal 		
MON								
TUE								
WED								
TUR								
FRI								
SAT								
SUN								
AVG								

My Self-Check Diary

	BREAKFAST				LUNCH			
WEEK	Pre-Meal 	Insulin Medication	Post-Meal 	Pre-Meal 	Insulin Medication	Post-Meal 		
MON								
TUE								
WED								
TUR								
FRI								
SAT								
SUN								
AVG								

ADVOCATE[®] *Redi-Code⁺*

Blood Glucose Meter

One Year Warranty Registration Card

Complete this warranty card and mail it promptly

Name _____

Street Address _____

APT# _____

City _____

State _____

ZIP _____

Country _____

Home Tel. _____

Meter Serial Number(See label on back of meter.)

Date of Purchase(M-D-Y) _____